

Accessories Ordering Information

| Item | Description | P/N |
|--------------------------------------|---|---------|
| BSQ-SP204S-A | An SP204S camera licensed for BeamSquared®. Sold as an accessory for those also purchasing a BSQ-XC130 | SP90631 |
| BSQ-Lens Kit UV-XNIR | Lens kit that includes 4 BeamSquared lenses: 500mm UV, 500mm VIS, 400mm NIR, 600mm XNIR | SP98009 |
| BSQ-Lens Kit 266-1550 | Lens kit that includes 5 BeamSquared lenses: 500mm UV, 500mm VIS, 400mm NIR, 400mm VIS, 400mm NIR, 400mm XNIR | SP90449 |
| BSQ-Lens Kit 650-1700 | Lens kit that includes 2 BeamSquared lenses: 400mm NIR, and 400mm XNIR. | SP90450 |
| BSQ-Lens UV 500mm | Single BeamSquared lens, 500mm focal length, A/R coated for 266-440nm | SP90451 |
| BSQ-Lens VIS 500mm | Single BeamSquared lens, 500mm focal length, A/R coated for 430-700nm | SP90452 |
| BSQ-Lens VIS 400mm | Single BeamSquared lens, 400mm focal length, A/R coated for 430-700nm | SP90453 |
| BSQ-Lens NIR 400mm | Single BeamSquared lens, 400mm focal length, A/R coated for 650-1000nm | SP90454 |
| BSQ-Lens XNIR 400mm | Single BeamSquared lens, 400mm focal length, A/R coated for 1000-1550nm | SP90455 |
| BSQ-Lens XNIR 600mm | Single BeamSquared lens, 600mm focal length, A/R coated for 1000-1550nm | SP90485 |
| BSQ-Lens UV 750mm | Single BeamSquared lens, 750mm focal length, A/R coated for 266-440nm | SP90554 |
| BSQ-Lens VIS 750mm | Single BeamSquared lens, 750mm focal length, A/R coated for 430-700nm | SP90555 |
| BSQ-Lens NIR 750mm | Single BeamSquared lens, 750mm focal length, A/R coated for 650-1000nm | SP90556 |
| BSQ-Lens XNIR 750mm | Single BeamSquared lens, 750mm focal length, A/R coated for 1000-1550nm | SP90557 |
| BSQ-Lens UV 1000mm | Single BeamSquared lens, 1000mm focal length, A/R coated for 266-440nm | SP90558 |
| BSQ-Lens VIS 1000mm | Single BeamSquared lens, 1000mm focal length, A/R coated for 430-700nm | SP90559 |
| BSQ SP300 or SP920 to SP204S Upgrade | Camera upgrade | SP90632 |
| BGS license for BSQ-SP204S | Includes BeamGage Standard software license in addition to BeamSquared software license | SP90633 |
| BGP license for BSQ-SP204S | Includes BeamGage Professional software license in addition to BeamSquared software license | SP90634 |
| BGP license for BSQ-XC130 | Includes BeamGage Professional software license in addition to BeamSquared software license | SP90508 |
| BSQ-XC130-KEY | Includes BeamSquared software license for XC-130 camera | SP90503 |
| BSQ - UV Reflective Filter | BeamSquared reflective ND2 filter, UV Grade Fused Silica, Inconel coating for 245-440nm | SP90568 |

